

XLR USB CHARGER



Shown with cell phone holder



Plugs conveniently into
the charger ports of
mobility scooters and
power chairs



Features and Benefits

- Charge your cell phone while on the go
- Can be used while power chair or scooter is in motion
- Charging detection—switches off automatically
- Standard with lanyard
- Short circuit, thermal and under-voltage protection
- XLR USB Charger will not fully discharge the battery

Compatible with: Most iPhone[®]*, Android and mobile devices

*iPhone[®] is a registered trademark of Apple, Inc.

5V DC, 1 A charger that is compatible with 24V and 36V systems