

# DELTA™ Ultra Light 1000, Full-Electric Low Bed

**New &  
Improved  
Game Changer**



15235



Hand crank located on the top of the footboard prevents bending for the caregiver when in the low position.




Truly universal with most head and footboards by other manufacturers. This makes managing inventory much easier!




Improved rail placement. Labels and color-coded springs assure accurate installation and patient safety.



Headboard taller than footboard. Better fits home decor!



Transition box is on the foot section NOT on the head and footboards. That's key to its interchangeability!



High-low shaft mounted on bed. Once installed, it can't be lost or forgotten on future installations!

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

**Ultra Light 1000 Full-Electric Low Bed**

15235	Ultra Light Full Electric Low Bed
15235BV-FR	Ultra Light Full Electric Low Bed with Full Length Side Rails (item #15001ABV)
15235BV-HR	Ultra Light Full Electric Low Bed with Half Rails (item #15201BV)
15235BV-PKG	Ultra Light Full Electric Low Bed with Full Length Side Rails (item # 15001ABV) and 80" Inner Spring Mattress (item #15006)
15235BV-PKG-T	Ultra Light Full Electric Low Bed with Full Length Side Rails (item #15001ABV) and 80" Therapeutic Support Mattress (item #15019)
15235BV-PKG-1	Ultra Light Full Electric Low Bed with Half Rails (item #15201BV) and 80" Inner Spring Mattress (item #15006)
15235BV-PKG-1-T	Ultra Light Full Electric Low Bed with Half Rails (item #15201BV) and 80" Therapeutic Support Mattress (item #15019)
15235BV-PKG-2	Ultra Light Full Electric Low Bed with Full Length Side Rails (item #15001ABV) and 80" Foam Mattress (item #15007)

Note: When the bed is placed in the "low" position, the side rails can not be lowered.

**FEATURES:**

- The Delta 1000 is truly universal because the headboards and footboards are completely interchangeable with Drive's and most other manufacturers bed ends, old and new.
- The transition box mounted on the foot section is key to interchangeability.
- The transition box allows you to change the rotation of the high-low shaft.
- Headboard is taller than footboard to better fit home decor.
- The wood grain panels are break and scratch resistant.
- Once installed, the high/low shaft will stay with the bed to eliminate loss or forgetting on future installations.
- Labels and color-coded springs assure accurate installation of side rails for patient safety.
- Unique motor is completely self-contained to reduce weight and noise.
- 9 1/2" to 23 1/2" deck height.
- Easy to set up.
- Top mounted Hi/Lo crank accessibility.
- 450 lb. weight capacity.
- Height adjustable motor comes pre-installed on the foot section.
- Caster easily moves from standard to low bed position.
- The motor assembly can be installed and removed with patient in bed.
- The motor is UL approved.
- In the event of a power failure, a 9 volt battery, located in the motor, can lower the head and foot sections nine times.
- New and improved hand pendant has large, easy to use controls.
- Bed ships in two cartons.
- Channel frame construction provides superior strength and reduced weight.
- Zinc coated spring deck.
- Extension kit available model #(15030EXTKIT) extends length from 80" to 84".
- Warranty: (Original purchaser)  
5 years on motor and hand control, Limited Lifetime on welds and frame,  
1 year on all other parts and components.

**SPECIFICATIONS**

Item #	Bed Height (Low Position) w/wheels	Bed Height (High Position) w/wheels	Overall Length	Overall Width	Product Weight	Carton Shipping Weight	Weight Limit
15035	9.5"	23.5"	88"	36"	167 lbs.	181 lbs.	450 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.